


HealingFit TWS

For detailed information on buttons and functions, please refer to the electronic manual available on the website or via the QR code

HealingFit TWS: Part Names and Functions



* Please check before use.

1. Charge before first use

* With the earbuds placed in the charging case, charge the case for at least 10 minutes.

2. Power on

* When you remove the earbuds from the charging case, they will turn on automatically (Voice prompt: "Power on")
* Alternatively, press and hold the touch button on the earbud for 5 seconds to turn it on.

3. Usage Modes

* The device supports both Dual Stereo Mode (using both earbuds) and Single Mode (using only one earbud).

4. HealingFit TWS Functions

HealingFit TWS provides three main functions: TES (Micro Electrical Stimulation) therapy, Sound Therapy, Bluetooth

* When using TES, wipe off any sweat or moisture from the electrode area and use the device only when it is completely dry.

* When charging with a charger other than a PC, please use a certified charger with a rated voltage of 5V.

* This device is not a medical device.

Pairing and Connection

1. Initial Pairing

* Entering pairing mode

When you remove both the left and right earbuds from the charging case, after about 3-5 seconds you will hear a "ding" sound and the earbuds will automatically enter pairing mode.

* Connecting on your phone

Go to Settings → Bluetooth on your phone, tap search, and select "HealingFit TWS" from the list. When you hear the voice prompt "Device connected," the connection is complete.

2. Reconnection and Connecting to a New Device

* Automatic reconnection: After the initial pairing, simply removing the earbuds from the charging case will automatically reconnect them to the previously paired smartphone.

* Important notes for connecting to a new device:

- To connect to a new device, make sure that HealingFit TWS is not currently connected to another device.
- If the earbuds automatically connect to a previous device, they may not appear in the search list on the new device. In this case, disconnect or turn off Bluetooth on the previous device to connect properly to the new one.

Earbud Usage Modes

* Dual Stereo Mode

- Features: Both the left and right earbuds are connected, allowing you to enjoy true stereo sound in both ears.
- Recommended: For music listening, it is recommended to use Dual Mode to experience the high-quality sound unique to Mobifren at its best.

* Single Mode

- Features: A mode that allows you to use only one earbud.
- Advantages: While one earbud is in use, the other can be charged in the charging case, enabling longer continuous use.

* Eartip Guide:

Three eartip sizes (Large · Medium · Small) are provided. Selecting the size that best fits your ears will help improve wearing comfort and sound quality.

Using Earphone Mode

Earphone Mode When the TES cable is disconnected

Music Modes

* HealingFit TWS offers five premium music modes, refined with MSTs tuning technology, which can be configured through the My HealingFit app.

* Choose the sound profile that suits your preferences to enjoy an immersive listening experience with optimal sound quality.

1. **MobiFren Mode** – Natural sound that stays true to the original audio

2. **BassBoost Mode** – Emotional sound with rich, enhanced bass

3. **Enhance Mode** – Clear sound with vivid detail

4. **Tender Mode** – Soft, comfortable sound that is gentle on the ears

5. **Space Mode** – Wide, immersive sound with a spacious, three-dimensional feel

* Music modes can be configured in the My HealingFit app

<p>1. Music Playback</p>	(In standby mode) L or R Double-tap quickly (2 times)	<p>2. Music Pause</p>	(While music is playing) L or R Quickly double-tap (2 times)	<p>3. Volume Up</p>	L or R Tap briefly at 1-second intervals Each tap increases the volume by 1 level
<p>4. Previous Track</p>	L button: Quickly triple-tap (3 times)	<p>5. Music Play</p>	R button Quickly triple-tap (3 times)	<p>6. Volume Down</p>	L or R Press and hold While holding, the volume decreases by 1 level at a time
<p>7. Voice Assistant (Sirj, etc.)</p>	(In standby mode) L or R Quickly tap 4 times	<p>8. ANC Function</p>	L or R After one short tap, press and hold Modes switch in the following order: → ANC ON → Transparency → ANC OFF	<p>ANC ON</p> Reduces external noise to block surrounding sounds. <p>Transparency</p> Captures external sounds and plays them along with music, allowing awareness of your surroundings. <p>ANC OFF</p> Used like regular earbuds with no noise control function.	
<p>9. Answer Call</p>	(When there is an incoming call) L or R Quickly double-tap (2 times)	<p>10. End Call</p>	(During a call) L or R Quickly double-tap (2 times)	<p>11. Reject Call</p>	(When there is an incoming call) L or R Press and hold for 2 seconds

Using Therapy Mode

You can use it more conveniently by installing the smartphone linkage app.

Scan the QR code or download from the App Store or Google Play
Search for "Healingfit"
Please install My HealingFit.

How to run the device after installing the APP

Make sure the device (mobile phone) on which the app is installed is connected to the Healing Fit TWS.

Move the screen according to the situation

- * When connected: Go to control page (play/pause therapy, change mode and adjust TES intensity)
- * When not connected: Go to the Bluetooth connection request page

[Troubleshooting]

If pairing fails, try restarting the device or resetting your phone's Bluetooth settings. If the app fails, restart the app or update to the latest version.

How to Use Therapy Mode (TES Therapy Mode)

* Starting Therapy Mode (Cable Connection)

- Preparation:** Keep the earbuds worn in both ears.
- Connection:** Open the charging case lid, take out the TES therapy cable, and connect it to both the left and right earbuds.
- Activation:** Once the cable is connected, Therapy Mode is automatically activated.

* Deactivating Therapy Mode and Duration

- Deactivation condition:** Therapy Mode will be deactivated 30 seconds after the cable is disconnected.
- Feature:** Even after the cable is removed, therapy continues for 30 seconds.

* Therapy Control Methods

- My HealingFit App:** Launch the app to conveniently adjust therapy intensity, volume, and duration.
- Earbud Touch Controls:** When the app is running, you can use the earbuds' touch controls to operate basic functions.

Therapy Mode (TES Cable Connected)

<p>1. Start Therapy</p>	(In standby mode) L or R Quickly double-tap (2 times)	<p>2. Pause Therapy</p>	(While music is playing) L or R Quickly double-tap (2 times)	<p>3. Volume Up</p>	L or R Tap briefly at 1-second intervals Each tap increases the volume by one level
<p>4. Decrease TES Intensity</p>	L or R: Tap once briefly, then press and hold While holding, the TES intensity decreases one level at a time	<p>5. Increase TES Intensity</p>	L or R: Tap once briefly, then press and hold While holding, the TES intensity increases one level at a time	<p>6. Volume Down</p>	L or R Press and hold While holding, the volume decreases by one level at a time

FAQ

Q. How should I adjust the intensity level of TES (microelectrical stimulation) to best use it?
A. Yes, this is normal. If stimulation occurs continuously in the same location, the contact resistance with the skin may change due to the skin condition of the contact area and the stimulation effect, and the intensity of the stimulation may vary. In this case, you can slightly adjust the position of the contact area (TES electrode) to the side. Try using it at a lower intensity level.

Q. What level of intensity is appropriate to use it?
A. Even with the same intensity, the degree to which it feels different varies from person to person, so it's impossible to make a general statement. We recommend using it at a lower intensity, starting from when you feel the stimulation. If you don't feel the stimulation with the electrical stimulation, you can use it while still feeling the stimulation. Even if you don't feel the stimulation, TES therapy still works. You can use it at whatever intensity you find comfortable.

Q. I set it to level 8 and am using it comfortably with a feeling of slight stimulation, but the stimulation suddenly started without me even changing the level. It feels hard. Is this normal?

A. Yes, this is normal. If stimulation occurs continuously in the same location, the contact resistance with the skin may change due to the skin condition of the contact area and the stimulation effect, and the intensity of the stimulation may vary. In this case, you can slightly adjust the position of the contact area (TES electrode) to the side. Try using it at a lower intensity level.

Product warranty

The product warranty period is one year from the date of purchase. Please refer to the website for details. Accessories are not repairable. If an accessory malfunctions within the warranty period, it will be serviced by the textbook and is not considered a reason for product exchange. Please back up important data in a safe place. We are not responsible for the loss of data stored on the product.

Product name	Model name		
Purchase date	Year	Month	Day
Where to buy			
Free warranty period:	Within 1 year from the date of purchase (excluding non-failure, marketing, and negligence) * Battery warranty period is 6 months.		

If a problem occurs while using the product, please check the user manual again and contact Mobifren Customer Service Center. We will promptly assist you.

AS Inquiries

For inquiries related to AS (after-sales service), as well as orders and delivery, please contact our Customer Service Center (1688-2251) for prompt assistance.

R-R-GTT-MFB-HC8000 / YU101532-21001B
Design & Technology by Mobifren KOREA, Assembled in China

For AS inquiries: 1688-2251

AS registration: www.mobifren.com/Customer support (39445)
88 Indong 52-gil, Gumi-si, Gyeongangbuk-do (Gupyeong-dong)
Mobifren Service Center TEL: 054-474-2220
FAX: 054-474-2251